

Lab Manual

NAME: _____

All rights reserved. © 2019 White Labs, Inc. All rights reserved. This book or any portion thereof may not be reproduced or used in any manner whatsoever without the express written permission of the publisher except for the use of brief quotations in a book review.



GENERAL GUIDELINES FOR ASEPTIC TECHNIQUE

Success microbiology labs rely heavily on good aseptic technique. Performing microbiological testing on dirty work surfaces or with nonsterile supplies can lead to issues. Aseptic technique is designed to provide a barrier between the environment and your culture. Below are some guidelines on performing good aseptic technique to reduce the potential of contaminating your cultures.

Clean Work Area


- Avoid air drafts and foot traffic
- Wash hands before touching sterile/sanitized equipment or containers
- Working surfaces should be washed with 70% isopropanol
- Be conscious of the work environment and equipment sterility/cleanliness

Good personal hygiene

- Tie back long hair
- Gloves should be worn and washed with 70% isopropanol
- Skin should be shielded with clothing or a lab coat

Handling

- Use sterilized glass or disposable plastic apparatuses when necessary
- Open and close caps on bottles, flasks, or tubes quickly to minimize air exposure
- Work with open bottles or flasks near a flame to avoid particles dropping into the apparatus




Serial Dilutions

MATERIALS

- Sterile water tubes (at least 2 will be needed for this dilution series)
- Transfer pipette
- Yeast slurry

METHOD

1. Set up test tubes in racks for dilutions with successive water tubes (depending on your desired dilution rate)
2. Label your tubes appropriately using a permanent marker.
3. Loosen the caps of the dilution tubes. Re-homogenize the yeast, carefully remove the lid, and pipette 1 mL of yeast and place into the 9mL water tube (for the primary dilution) using sterile transfer pipette. Pump the pipette to mix up the yeast and water. You have just made a 1:10 dilution
4. Using the same pipette, remove 1 mL of diluted yeast and place into the next water tube. This will make a 1:100 dilution.



HLP TESTING

MATERIALS

- Diluted yeast slurry or beer sample
- 15 mL conical tube
- Prepared HLP media (Microwaved)
- 70% Isopropanol
- Transfer pipette
- Bunsen Burner

METHOD

1. Wipe down working space with 70% isopropanol and turn on Bunsen burner
2. Label 15ml conical tubes with your sample name
3. Arrange tubes as close as possible next to flame to work inside the “aseptic bubble” generated by the flame.
4. Aseptically transfer 1ml of sample from to the sterile conical tube. Replace cap after dispensing sample.
5. Fill tube with approximately 14 ml of HLP media or until the media reaches the top of the tube. After securely closing the cap, invert tube multiple times to thoroughly mix sample.
6. After 5 days inspect plates. Negative tubes have no growth. Positive samples have white dots or tubes are opaque. (See included positive control)

STREAK PLATE


MATERIALS

- Agar plate with colony that you'd like to streak out (working plate)
- Sterile one time use inoculating loop(s)
- Pre-poured agar plate or slant (Generally YPD will be used for *Saccharomyces cerevisiae*)

METHOD

*Note: There are several acceptable methods that can be utilized to streak plates. One such method is described below.

1. Wipe down working space with 70% isopropanol and turn on Bunsen burner
2. Label YPD plate with your name, group number, and date
3. Arrange plates as close as possible next to flame to work inside the "aseptic bubble" generated by the flame.
4. Pick up 1 colony from working plate
5. Streak in a zigzag on 1/3 of the plate
6. Use new loop to streak through first zigzag
7. Use new loop to streak through second zigzag
8. Incubate at 30°C for 3-5 days



FORCED WORT

MATERIALS

- Wort (with out yeast)
- Yeast slurry
- Sterile container- Whirlpak bag, sanitized flask or mason jar with foil
- Warm incubation space

METHOD

1. Sanitize sample valve
2. Aseptically transfer 50mL of wort to sterile container. Run wort out of the sample valve for at least 5 seconds before collecting a sample.
 - a. Option to pipette sample from fermentor
3. Close container and incubate in warm area

Results:

Clear = beer is clean

Cloudy wort or wort with bubbles = contamination

Duration	Result
1 day	Very dirty, clean heat exchanger and hoses. Beer will need to be dumped.
2-3 days	Major contamination. Need to clean problem, beer most likely will be affected. Do not collect yeast for re-use from this batch.
3-6 days	Mild contamination build up, clean problem. Beer may or may not be affected.
7 or more	Very clean, keep up the good work

FORCED DIACETYL

MATERIALS

- Two glasses
- Aluminum Foil
- Hot Water Bath
- Thermometer

METHOD

Procedure:

1. Heat the water bath to 140 to 160 F
2. Collect beer into each glass and cover with aluminum foil
3. Place one glass in the hot water bath, while keeping the other at room temperature. Do not microwave!
4. After 10 to 20 min remove the beer from the hot bath, and cool to the same temperature as the other sample. Ice water bath is effective for cooling.
5. Remove the aluminum foil and smell each sample. If you smell the buttery character of diacetyl in either or both samples, you know your beer has the diacetyl precursor.

CELL COUNTING & VIABILITY


MATERIALS

- Methylene blue solution
- Hemocytometer
- Hand held counter
- Transfer pipets
- Kimwipes
- Gloves (optional)

You will need a microscope with 400X capability (40x lens + 10x eyepiece).


METHOD

1. Make certain hemocytometer is clean and dry before use. Hemocytometers can be easily cleaned with water and/or 70% EtOH. Chamber and cover slip may be scrubbed gently using a lint free towelette (kimwipe).
2. Dilute yeast sample to an appropriate concentration. Usually a 1:100 or 1:1000 works well for yeast slurry. Sample can be diluted with deionized/distilled water or with 0.5% H₂SO₄ if cells clump excessively. **Note your dilution factor.**
3. Invert sample several times to mix, taking care not to introduce bubbles.
4. Mix 1mL of your diluted yeast sample with 1mL of methylene blue solution and allow to incubate for 1 - 2 minutes.
5. Invert sample 1-2 times and take up sample by using a transfer pipette.
6. Fill chamber by placing 2 drops of sample onto center of counting grid.
7. Gently position cover slip so that glass covers both counting areas equally.
8. Carefully place hemocytometer on microscope stage. As you focus at each objective leading up to the 40X lens, note the distribution of yeast cells. If cells are well distributed, then you can use the short cell count method.



CONSIDERATIONS FOR WHEN CELL COUNTING


1. You will be counting squares within the 1mm² ruled area centrally located on the chamber (see images below, in next section)
2. It is helpful to establish a counting protocol for all cell counts and stick with it! Consistency is key!
--For example, cells touching or lying on the top and right boundary lines are not counted, whereas cells touching or lying on the bottom or left boundary lines are counted (see image below).
3. Yeast cell buds emerging from mother cells are counted as a separate cell if the bud is at least one-half the size of the mother cell.
4. If performing viability counts, dead cells will stain dark blue. Non-viable cells do not have the metabolic capability to reduce the blue dye.
5. Do not count cells that are pale blue in color as dead. Some budding cells will also stain blue, do not count these cells as non-viable.
6. If you are performing a cell count and viability count simultaneously, it is best to count all cells on the hand-held counter and record noted dead cells on a written tally.



SHORT METHOD: FOR AN EVENLY DISTRIBUTED SAMPLE

1. You will be counting cells within the 5 numbered squares (see image below)
2. Estimate the total number of cells in grid by taking # of cells counted in 5 numbered squares and multiplying by 5
3. Your yeast count equation is:

$$\text{Yeast cells/mL} = \text{Total number of cells in grid} \times \text{dilution factor} \times 10^4 \text{ (or } 10,000, \text{ this number is a constant.)}$$



Example: Let's say you counted 136 cells in the 5 numbered squares, using a 1:200 dilution.

Your cell count would be: $136 \times 5 \times 200 \times 10,000 = 1360000000 = 1.36 \times 10^9 = 1.36 \text{ billion cells/mL}$

YEAST PITCHING CALCULATION

Example:

Parameter	Value	Unit
Batch size	5	Gallon
Strength of wort	1.056	S.G
Wanted re-pitching rate	1	Million cells / ml / degree Plato
Concentration of slurry	0.8	Billion cells / ml

Question: How much yeast slurry (volume) should the brewer add to achieve the wanted pitching rate?

First, let's calculate the total number of cells we want to add to the full batch of wort

Goal: We want 1 million cells / ml / degree Plato. With a 1.056 S.G wort, we need to convert to Plato. (Use chart) Therefore we want 13.8 million cells / ml or $13.8 * 10^6$ cells/ml

We have 5 gal wort. 5 gal is 18927mls.

We have just calculated that we have 18927 ml wort. Therefore, our total goal of how many cells to add is:

$$18927ml * 13,800,000 \text{ cells/ml} = 2.6 * 10^{11} \text{ cells}$$

Slurry concentration obtained by cell count is:

$$800 \text{ million cells / ml, or } 8 * 10^8 \text{ cells /ml}$$

How much slurry to add?

Our goal was to add $2.6 * 10^{11}$ cells to the fermentor and we have 0.8 million cells/ml in our slurry. To calculate this we must use the formula for pitching

$$\text{Volume of slurry to pitch (ml)} = \frac{\text{total cells wanted (cells)}}{\text{concentrations of slurry (cells/ml)}}$$

So, now we have all we need:

$$\text{Volume to pitch (ml)} = \frac{2.6 * 10^{11} \text{ cells}}{8 * 10^8 \text{ cells/ml}}$$

$$\text{Volume to pitch (ml)} = 325\text{ml}$$

SPECIFIC GRAVITY TO PLATO CONVERSION FIGURE


Brix	Plato	SG	Brix	Plato	SG	Brix	Plato	SG
0.0	0.0000	1.0000	13.4	13.4027	1.0543	26.8	26.7948	1.1140
0.2	0.1970	1.0008	13.6	13.6028	1.0551	27.0	26.9944	1.1150
0.4	0.3970	1.0016	13.8	13.8029	1.0560	27.2	27.1940	1.1159
0.6	0.5970	1.0024	14.0	14.0030	1.0568	27.4	27.3936	1.1168
0.8	0.7970	1.0031	14.2	14.2030	1.0577	27.6	27.5932	1.1178
1.0	0.9970	1.0039	14.4	14.4031	1.0586	27.8	27.7928	1.1187
1.2	1.1970	1.0047	14.6	14.6031	1.0594	28.0	27.9924	1.1197
1.4	1.3971	1.0054	14.8	14.8032	1.0603	28.2	28.1919	1.1206
1.6	1.5971	1.0062	15.0	15.0032	1.0611	28.4	28.3915	1.1216
1.8	1.7971	1.0070	15.2	15.2033	1.0620	28.6	28.5910	1.1225
2.0	1.9972	1.0078	15.4	15.4033	1.0628	28.8	28.7905	1.1235
2.2	2.1972	1.0086	15.6	15.6033	1.0637	29.0	28.9901	1.1244
2.4	2.3973	1.0094	15.8	15.8034	1.0646	29.2	29.1896	1.1254
2.6	2.5973	1.0101	16.0	16.0034	1.0654	29.4	29.3891	1.1263
2.8	2.7974	1.0109	16.2	16.2034	1.0663	29.6	29.5886	1.1273
3.0	2.9975	1.0117	16.4	16.4034	1.0672	29.8	29.7880	1.1282
3.2	3.1975	1.0125	16.6	16.6034	1.0680	30.0	29.9875	1.1292
3.4	3.3976	1.0133	16.8	16.8034	1.0689	30.2	30.1870	1.1302
3.6	3.5977	1.0141	17.0	17.0034	1.0698	30.4	30.3864	1.1311
3.8	3.7977	1.0149	17.2	17.2034	1.0706	30.6	30.5859	1.1321
4.0	3.9978	1.0157	17.4	17.4034	1.0715	30.8	30.7853	1.1330
4.2	4.1979	1.0165	17.6	17.6034	1.0724	31.0	30.9847	1.1340
4.4	4.3980	1.0173	17.8	17.8034	1.0733	31.2	31.1841	1.1350
4.6	4.5981	1.0181	18.0	18.0033	1.0741	31.4	31.3835	1.1359
4.8	4.7982	1.0189	18.2	18.2033	1.0750	31.6	31.5829	1.1369
5.0	4.9983	1.0197	18.4	18.4033	1.0759	31.8	31.7823	1.1379
5.2	5.1984	1.0205	18.6	18.6032	1.0768	32.0	31.9817	1.1389
5.4	5.3985	1.0213	18.8	18.8032	1.0777	32.2	32.1810	1.1398
5.6	5.5986	1.0221	19.0	19.0031	1.0785	32.4	32.3804	1.1408
5.8	5.7987	1.0229	19.2	19.2030	1.0794	32.6	32.5797	1.1418
6.0	5.9988	1.0237	19.4	19.4030	1.0803	32.8	32.7791	1.1428
6.2	6.1989	1.0245	19.6	19.6029	1.0812	33.0	32.9784	1.1437
6.4	6.3990	1.0253	19.8	19.8028	1.0821	33.2	33.1777	1.1447
6.6	6.5991	1.0261	20.0	20.0027	1.0830	33.4	33.3770	1.1457
6.8	6.7992	1.0269	20.2	20.2026	1.0839	33.6	33.5763	1.1467
7.0	6.9994	1.0277	20.4	20.4025	1.0848	33.8	33.7756	1.1477
7.2	7.1995	1.0285	20.6	20.6024	1.0857	34.0	33.9749	1.1487
7.4	7.3996	1.0294	20.8	20.8023	1.0866	34.2	34.1741	1.1497
7.6	7.5997	1.0302	21.0	21.0021	1.0875	34.4	34.3734	1.1507
7.8	7.7998	1.0310	21.2	21.2020	1.0884	34.6	34.5727	1.1516
8.0	7.9999	1.0318	21.4	21.4018	1.0892	34.8	34.7719	1.1526
8.2	8.2000	1.0326	21.6	21.6017	1.0901	35.0	34.9711	1.1536
8.4	8.4002	1.0334	21.8	21.8015	1.0911	35.2	35.1703	1.1546
8.6	8.6003	1.0343	22.0	22.0014	1.0920	35.4	35.3695	1.1556
8.8	8.8004	1.0351	22.2	22.2012	1.0929	35.6	35.5687	1.1566
9.0	9.0005	1.0359	22.4	22.4010	1.0938	35.8	35.7679	1.1576
9.2	9.2006	1.0367	22.6	22.6008	1.0947	36.0	35.9671	1.1586
9.4	9.4007	1.0376	22.8	22.8006	1.0956	36.2	36.1663	1.1596
9.6	9.6009	1.0384	23.0	23.0004	1.0965	36.4	36.3655	1.1606
9.8	9.801	1.0392	23.2	23.2002	1.0974	36.6	36.5646	1.1617
10.0	10.0011	1.0400	23.4	23.4000	1.0983	36.8	36.7638	1.1627
10.2	10.2012	1.0409	23.6	23.5997	1.0992	37.0	36.9629	1.1637
10.4	10.4013	1.0417	23.8	23.7995	1.1001	37.2	37.1620	1.1647
10.6	10.6014	1.0425	24.0	23.9992	1.1011	37.4	37.3612	1.1657
10.8	10.8015	1.0434	24.2	24.1990	1.1020	37.6	37.5603	1.1667
11.0	11.0016	1.0442	24.4	24.3987	1.1029	37.8	37.7594	1.1677
11.2	11.2017	1.0450	24.6	24.5984	1.1038	38.0	37.9585	1.1688
11.4	11.4018	1.0459	24.8	24.7982	1.1047	38.2	38.1576	1.1698
11.6	11.6019	1.0467	25.0	24.9979	1.1057	38.4	38.3566	1.1708
11.8	11.8020	1.0475	25.2	25.1976	1.1066	38.6	38.5557	1.1718
12.0	12.0021	1.0484	25.4	25.3972	1.1075	38.8	38.7548	1.1728
12.2	12.2022	1.0492	25.6	25.5969	1.1084	39.0	38.9538	1.1739
12.4	12.4023	1.0501	25.8	25.7966	1.1094	39.2	39.1529	1.1749
12.6	12.6024	1.0509	26.0	25.9963	1.1103	39.4	39.3519	1.1759
12.8	12.8025	1.0518	26.2	26.1959	1.1112	39.6	39.5509	1.1770
13.0	13.0026	1.0526	26.4	26.3956	1.1122	39.8	39.7500	1.1780
13.2	13.2027	1.0534	26.6	26.5952	1.1131	40.0	39.9490	1.1790



EXERCISES



CELL COUNTING EXERCISE

Counting the amount of cells dead and alive in each photo or microscope view.

This view is at 100X - You can see all of the 25 squares on the hemocytometer







Based on the following parameters, calculate the cell count of the slurry.

Dilution Factor 1: 200

Yeast cells/mL = Total number of cells in grids (reminder to multiple by 5 for the average) x dilution factor x 10^4 (or 10,000, this number is a constant.)

PITCHING RATE EXERCISE

Parameter	Value	Unit
Batch size	5	Gallon
Strength of wort	1.050	S.G
Wanted re-pitching rate	1	Million cells / ml / degree Plato
Concentration of slurry	Use concentration from cell counting exercise	Billion cells / ml

Question: How much yeast slurry (volume) should the brewer add to achieve the wanted pitching rate? Reminder of the equation is below

$$\text{Volume of slurry to pitch (ml)} = \frac{\text{total cells wanted (cells)}}{\text{concentrations of slurry (cells/ml)}}$$

